

AGEX^{fill}

FEATURES

- REGENERATION
- ELASTICITY
- TISSUE REDENSIFICATION
- SOFTNESS AND FLEXIBILITY
- RECEPTIVE STIMULATION

Filler that has a biostimulating action allowing the complete rejuvenation of the skin. Restores elasticity, hydration and stimulates cell renewal.

Included in the packaging: 27 G x 37 mm cannula
26 G x 13 mm needle

Also compatible with:
27 G x 13 or 19 mm needles, up to 30 G x 13 mm

PROTOCOL	2-3 treatments every 2-3 months
DURATION	2-3 months
COMPOSITION	Crosslinked hyaluronic acid + Linear hyaluronic acid
DOSAGE	1.1 ml
MOLECULAR WEIGHT	2 M Da + 1 M Da Crosslinked HA + 10% Linear HA 500,000 Da
INJECTION POINT	Dermi
CONCENTRATION	25 mg/ml
PACKAGING	3 pre-filled single-use syringe
CROSS-LINKING AGENT	Residual BDDE < 0,1 ppm
ELASTIC MODULUS G'	30 Pa
VISCOUS MODULUS G''	10 Pa

CROSSLINKING



VOLUME



BIOSTIMULATION



- FACE
- NECK/ DÉCOLLETÉ
- BACK OF THE HANDS



VIDEO
ON FILLER
PROCEDURE



AGEX^{fill} VOLUME

FEATURES

- VOLUME
- PLUMPNESS
- RESTRUCTURING

Indicated to correct the average volumes of the face. It adapts to the individual morphology, respecting the naturalness of the traits.

Included in the packaging: 27 G x 13 mm needle
27 G x 19 mm needle

Also compatible with:
25 G x 38 mm / 50 mm cannula

DURATION	5/6 months
COMPOSITION	Crosslinked hyaluronic acid + Linear hyaluronic acid
DOSAGE	1 ml
MOLECULAR WEIGHT	2 M Da + 1 M Da Crosslinked HA + 8% Linear HA 1 M Da
INJECTION POINT	Hypodermis
CONCENTRATION	25 mg/ml
PACKAGING	1 pre-filled single-use syringe
CROSS-LINKING AGENT	Residual BDDE < 0,1 ppm
ELASTIC MODULUS G'	100 Pa
VISCOUS MODULUS G''	20 Pa

CROSSLINKING



VOLUME



BIOSTIMULATION



- MEDIUM VOLUME RESTORATION
- CHEEKS
- ASYMMETRY CORRECTION
- MALAR REGION
- NASOLABIAL WRINKLES
- MEDIUM DEPTH MARIONETTE WRINKLES
- VOLUMETRIC LIP AUGMENTATION



VIDEO
ON FILLER
PROCEDURE



AGEX^{fill} ULTRA

FEATURES

- RESTORATION OF DEEP VOLUMES
- CORRECTION OF DEEP SKIN SAGGING

Ideal for restoring the deep volumes of the face and correcting deep skin sags.

Included in the packaging: 25 G x 38 mm cannula
23 G x 19 mm needle

Also compatible with: 27 G x 13 mm needle

<u>DURATION</u>	5/6 months
<u>COMPOSITION</u>	Crosslinked hyaluronic acid + Linear hyaluronic acid
<u>DOSAGE</u>	1 ml
<u>MOLECULAR WEIGHT</u>	2 M Da + 1 M Da Crosslinked HA + 5% Linear HA 500.000 Da
<u>INJECTION POINT</u>	Hypodermis
<u>CONCENTRATION</u>	25 mg/ml
<u>PACKAGING</u>	1 pre-filled single-use syringe
<u>CROSS-LINKING AGENT</u>	Residual BDDE < 0,1 ppm
<u>ELASTIC MODULUS G'</u>	200 Pa
<u>VISCOUS MODULUS G''</u>	40 Pa



- RINOFILLER
- SEVERE NASOLABIAL FOLDS
- SEVERE MARIONETTE WRINKLES
- REDEFINITION OF MANDIBULAR CONTOURS
- MALAR REGION
- FRONTAL REGION
- ZYGOMATIC REGION
- TEMPORAL REGION
- MEDICAL MENTOPLASTY

CROSSLINKING



VOLUME



BIOSTIMULATION



VIDEO
ON FILLER
PROCEDURE

